

# Gymna 400 series

---

Feel the Power



# Gymna 400 series

## Feel the power

### How it works


*Treatment with electrotherapy is well accepted in the field of physical therapy and it involves the introduction of some physical energy into a biologic system. This energy brings about one or more physiological changes that are used for therapeutic benefit.*



### Benefits

- Advanced healing process
- More effective circulatory system
- Increased muscle tone
- Increased speed and strength of muscle contractions
- Increased rate of absorption
- Reduction in pain


#### Evidence of TENS for pain management






# Gymna 400 series

## 3 powerful goals




PAIN RELIEF

*Efficient &  
Focused*



CELLULAR EFFECTS

*Effective &  
Flexible*



MUSCLE STIMULATION

*Dynamic &  
Fast*




# Gymna 400 series

## Flexibility to fit your needs

### Colour-guided therapy

Every therapy has its own colour for recognition at a glance.



### Electrode placement with real-life pictures

Giving guidance on positioning electrodes, vaco cups, laser or ultrasound heads.



### Customisation wizard

Set the preferred user interface and customise the device for your practice.



### Full screen parameter view

Therapy screen with possibility to show parameters full screen for progress view at a glance.





# Gymna 400 series

## Our expertise your gain

### *Guided Therapy system*

The building blocks for GTS are our more than 40 years experience in electrotherapy, the expertise we have built up over these years and our constant push to improve the ease of use.

GTS helps and guides the beginning therapist in choosing the best treatment for each individual patient and supports the experienced therapist with its ease of use and flexibility.

#### **Main GTS features:**

- Intuitive navigation
  - Colour-guided therapy*
  - Uniformity leads to fast action & response*
- Medical E-book
  - Anatomical library*
  - Clinical information screens*
  - Help screens*
- Evidence-based therapy
  - In line with latest scientific studies*
  - New protocols can be uploaded via USB*



# Gymna 400 series

## Elegance and ease of use combined

### *The new Gymna 400 ViP line*

*The new Gymna 400 ViP offers a total and elegant solution for your practice, consisting of a Gymna 400 electrotherapy device, a vaco device and the new Gymna Mobile 400.*

*The new Gymna Mobile 400 is our latest solution to help you in your everyday practice. A clever mobile storage system for accessories and cables with immediate and easy access when necessary.*





# Gymna 400 series

## Gymna Mobile 400



*Laser probes*



*Vacuum cups,  
holder & cables*



*Electrodes &  
Electrode cables*




*Tissue box*



*Ultrasound  
heads*



*Vaco 400 device*



# Gymna 400 series

## Product offering

### Combi 400 ViP

Combination device with electrotherapy, ultrasound, laser (optional), including modular vaco unit and Gymna Mobile 400.



### Combi 400V

Combination device with electrotherapy, ultrasound, laser (optional), including modular vaco unit.



### Combi 400M

Combination device with electrotherapy, ultrasound, laser (optional), including Gymna Mobile 400.



### Combi 400

Combination device with electrotherapy, ultrasound, laser (optional).



### Pulson 400M

Device for ultrasound applications. Including 2 multi-frequency ultrasound heads (1 and 3 MHz), 1 and 4 cm<sup>2</sup> and Gymna Mobile 400.



### Pulson 400

Device for ultrasound applications. Including 2 multi-frequency ultrasound heads (1 and 3 MHz), 1 and 4 cm<sup>2</sup>



### Duo 400 ViP

Device for electrotherapy applications, including modular vaco unit and Gymna Mobile 400.



### Duo 400V

Device for electrotherapy applications, including modular vaco unit.



### Duo 400M

Device for electrotherapy applications, including Gymna Mobile 400.



### Duo 400

Device for electrotherapy applications.



### Vaco 400

Modular vacuum unit with 2 independent channels, controlled via Combi 400 or Duo 400



### Mobile 400

Mobile storage system for accessories and cables with immediate and easy access when necessary.



# Gymna 400 series

## Accessories 400 series

Electrode cables (2 pcs)

Rubber electrodes (2 pcs) 2x

Chamex bags 6x8 cm (4 pcs)

Fixing straps 5 x 60 cm (4x)

US treatment head 4 cm<sup>2</sup> (1 & 3 MHz)

US treatment head 1 cm<sup>2</sup> (1 & 3 MHz)

VAS-score

Touch screen pen

Contactgel 500ml

Test plug

Vacuum tube light grey (2 pcs)

Vacuum tube dark grey (2 pcs)

Vacuum electrode 60 mm (4 pcs)

Sponge for vacuum electrodes (4 pcs)

Monoprobe 400

Clusterprobe 400

Vaco 400

Gymna Mobile 400

	COMBI				DUO				PULSON	
	400VIP	400V	400M	400	400VIP	400V	400M	400	400	400M
	■	■	■	■	■	■	■	■		
	■	■	■	■	■	■	■	■		
	■	■	■	■	■	■	■	■		
	■	■	■	■	■	■	■	■		
	■	■	■	■					■	■
	○	○	○	○					■	■
	■	■	■	■	■	■	■	■	■	■
	■	■	■	■	■	■	■	■	■	■
	■	■	■	■					■	■
	■	■	■	■	■	■	■	■		
	■	■			■	■				
	■	■			■	■				
	■	■			■	■				
	■	■			■	■				
	○	○	○	○						
	○	○	○	○						
	■	■	○	○	■	■	○	○		
	■	○	■	○	■	○	■	○	○	■

■ Standard

○ Optional

# Gymna 400 series

## Characteristics

### Direct therapy modes:

Electro current forms	31
Iontophoresis	18
Ultrasound & Phonophoresis	2
Combination therapy	18
Laser therapy	2

### Evidence based protocols:

Objectives	190
Indications	279
Body Area	234
Cellular effects	16
Diagnostics	12

### Memory free locations:

Favourites & own programs	500
Diagnostics	200
Own sequential	100
Shared programs	50





# Gymna 400 series

## Technical Data 400 series

[illegible]

## GYMNA 400 SERIES

*Choose your  
preferred colour*



Black



White

**BUTTONS**  
*available  
in two colours*



Silver

Black

### Cryotherapy



### Tecar - Diathermy



### Shockwave Therapy



### Physio Care



### Treatment Tables



### Electrotherapy



### Active Motion



Your dealer

[gymna.com](http://gymna.com)

**gymna**<sup>®</sup>

GymnaUniphy NV Pasweg 6A | B-3740 Bilzen, Belgium

Tel: +32 (0)89 510 510 | Fax: +32 (0)89 510 511